

RS Leucocyte Removal Filter



For Salvaged Blood

- Clinically proven media technology
- Reduces leucocytes, fat and microaggregate debris found in salvaged blood
- Designed for use with either washed or unwashed intra-operative salvaged blood
- Low filter hold up (31 mL after recovery) and vented spike allowing excellent blood product recovery without the need for saline flushing
- High flow rate for rapid infusion
- Automatic priming filter and self-levelling administration set designed for ease of use
- Latex free

Performance Characteristics

Leucocytes and lipid particles were studied from both washed and unwashed salvaged blood, obtained from patients undergoing coronary artery bypass procedures and valve replacements or repairs.

Pre- versus post-filtration values were analysed using a paired t-test. The results demonstrate significant reductions of both leucocytes (99%; $P < 0.001$) and lipid particles (82%; $P < 0.001$) from salvaged blood.

| | Average Volume (range) | HCT % (range) | Flow Rate mL/min (range) | Leucocytes/mL | | Lipid Particles/mL* | |
|------------------------|---------------------------|---------------|--------------------------|----------------------------|-----------------|---------------------|-----------------|
| | | | | Pre-Filtration | Post-Filtration | Pre-Filtration | Post-Filtration |
| Unwashed Blood (n = 8) | 678 mL (450 – 1000 mL) | 20 (16 – 23) | 52 (3 – 200) | 18,385 (5,280 – 67,200) | 50 (0.1 – 224) | – | – |
| Washed Blood (n = 15) | 1 bowl (1 – 2 bowls) | 61 (53 – 68) | 82 (41 – 112) | 23,233 (8,200 – 39,200) | 189 (0.5 – 750) | 260 (12 – 869) | 48 (12 – 167) |

*Lipid particles < 10 μ m. There were no detectable lipid particles > 10 μ m.

Lipid was calculated as the product of fat concentration (globules/mL) and volume of blood product (mL).

Ordering Information

| Reorder Code | Description | Pkg |
|--------------|---|---------|
| RS1VAE | RS Leucocyte Removal Filter for Salvaged Blood with Self Levelling Administration Set | 20/case |



Information contained herein for distribution outside of the USA only.

© 2012, 2013 Haemonetics Corporation. Haemonetics and Haemonetics The Blood Management Company are trademarks or registered trademarks of Haemonetics Corporation in the USA, other countries, or both. All rights reserved. 04.2013 USA. COL-COPY-000576-IE(AA)

Haemonetics S.A. ■ Signy Centre ■ Rue des Fléchères ■ P.O. Box 262 ■ 1274 Signy-Centre ■ Switzerland ■ www.haemonetics.com