MCS®+ 9000

The Added Advantage In Platelet Collection





Add More Platelets in More Places with the MCS®+

Donors

Collecting more platelets starts with reaching more donors. And many donors—from smaller individuals to those with lower platelet counts or time constraints—are more suited to single unit collections. Many other potential donors only consider donation when mobile units come to them. So optimize all collections with the proven technology and unique mobility of the MCS+ device for every donor everywhere.

Places

- Fixed Sites
 - Collect single platelet units more cost-effectively
 - Increase comfort for all donors using saline return
- Mobile Sites
 - Access new platelet donors and build your donor base
 - Increase mobile drive profits by collecting more platelets



Platelets

Mobile	Collection	System
Only 56 pounds for ultimate portability—lightest in the industry	Cost-effective collection from single platelet donors	Automated pressure cuff, saline return and reduced citrate usage for donor comfort ^{1,2}
 Rugged device withstands any mobile environment 	Proven separation technology State-of-the-art leukoreduction	Simple installation, operation and automated cycles for operator ease of use
 Customized transport case doubles as a device stand 	with continuous filtration Double unit and concurrent plasma capabilities	Robust family of devices with quiet performance

Partnering with Haemonetics

- Industry-leading equipment reliability and disposable quality 300%+ better than the general average,³ optimizing uptime and increasing efficiency
- Expert support for program implementation and extensive post-installation monitoring delivers enhanced staff productivity to meet collection and operational goals

¹ Only platelet device with zero vasovagal or citrate reactions per 1,000 donors. M Kim, Korean Red Cross, Seoul; J Shin, Soonchunhyang University Hospital, Seoul; et al. A Comparative Study Among Three Cell Separators. Transfusion 2002; 42:35S.

^{2 &}quot;Donors sensitive to citrate could benefit from the reduced citrate infusion rates of MCS+." S M Radojska, S M Picker, et al, Medizin, University of Cologne, Germany. Prospective Comparison of Citrate Donor Load Between Different Multicomponent Apheresis Systems. Transfusion 2005; 45:47A.

 $^{3~}M~Harry, R~Schroeder.~1999, \\ "Six~Sigma-Breakthrough~Management~Strategy."$

Bring Safety, Efficiency and Cost-Effectiveness to Platelet Collection

Why Single Donor Platelet Collection?

Clinical Advantages

In comparison with whole blood random donor platelets (RDP), single donor platelets (SDP) offer several significant advantages:

- Provides a consistent dose for transfusion.
- Delivers less risk of bacterial contamination⁴
- Exhibits less cellular activation⁵
- Ensures better platelet recovery rates and longer survival⁶
- Decreases risk of infections, allergic reactions and platelet alloimmunization for patients

Processing Advantages

- SDP collection requires only one donor per dose, while RDP collection requires multiple donors.
- SDP collection eliminates the need for additional processing and testing, providing platelets that are ready faster to meet patient needs.
- Software adjusts the plasma levels each cycle, as needed, to achieve a consistent target yield and concentration.



4 Ness et al. Transfusion 2001;41:857-861. Yomtovian et al. Transfusion 2006;46:719-730. 5 M Bock et al. Transfusion Medicine 2002;12:317-324. 6 D Arnold. N Heddle. Transfusion 2006:46:257-264.



Studies Indicate the Risk of Bacterial Contamination Is Up to 5x Less with SDP					
	RDP	SDP	Р		
Ness et al.	1:4818	1:15,078	N/A		
Yomtovian et al.	1:418 (pools)	1:2213	0.001		
Yomtovian et al. (contaminated units transfused)	1:294	1:8011	0.001		
Ness et al. Transfusion 2001;41:857-861. Yomtovian et al. Transfusion 2006;46:719-730.					

Apheresis SDP Had 19% Better Recovery and 33% Longer Survival than WB RDP					
In Vitro 5-day	Apheresis	RDP	Р		
P-Selectin % Morphology Score	16.4 257	22.0 232	0.035 0.008		
In Vivo	Comparison of recovery and survival of whole blood platelets vs. apheresis platelets				
	Whole Blood	Apheresis	Р		
% Recovery Survival (hrs)	43 101	51 134	0.03 <0.001		
D Arnold, N Heddle. Transfusion 2006;46:257-264.					

To find out more about how Haemonetics continues to advance healthcare with innovative blood management solutions, call us at 1-800-225-5242.

GLOBAL HEADQUARTERS

Haemonetics Corporation 400 Wood Road Braintree, Massachusetts 02184 USA

www.haemonetics.com